

**Bioactive Glass creates  
the state of healing in the root canal.**



Root Canal Sealer  
**NISHIKA CANAL SEALER BG**  
/// Bioactive Glass

 **NISHIKA** Nippon Shika Yakuhin Co.,Ltd.

URL: <http://www.nishika.co.jp/english/>

*Bioactive Glass induces the optimal condition for periapical healing.*

**Next-generation Root Canal Sealer  
“NISHIKA CANAL SEALER BG”**



**Chiaki Kitamura**  
DDS, PhD  
Professor and Chairperson  
Div. of Endodontics and Restorative Dentistry,  
Dept. of Oral Functions, School of Dentistry  
Kyushu Dental University



**Ayako Washio**  
DDS, PhD  
Assistant Professor,  
Div. of Endodontics and Restorative Dentistry,  
Dept. of Oral Functions, School of Dentistry  
Kyushu Dental University

**Introduction**

Three requirements for root canal sealer	02
Application of bioceramics to root canal sealer	02
Characteristics of Bioactive Glass	02
New Bioceramic root canal sealer	03

**Features of NISHIKA CANAL SEALER BG**

Physicochemical properties	04
Excellent biocompatibility	04
Excellent sealing ability	05
Handleability	06
Removability	06
Precaution	07

**References**

References	07
------------	----

**Presentation**

Packaging	07
Composition	07

**Introduction****Three requirements for root canal sealer**

Conventional root canal sealers as typified by Grossman's formula were specifically for filling the space between gutta-percha points and root canal walls, but hardly ideal since neither adhesive nor bonding effect occurs.

Recently, a root canal sealer that strongly bonds to root canal walls and exerts high sealing properties has been desired. Root canal sealer is also required to have higher biocompatibility so as not to induce inflammatory reactions in periapical tissue. Furthermore, usability is a significant factor in clinical applications. In short, root canal sealer is required to have the three properties of excellent biocompatibility, sealing ability and removability.

**Application of bioceramics to root canal sealer**

Mineral Trioxide Aggregate (MTA) offers nearly ideal performance as a root canal filling material. MTA has triggered the development of new bioceramic root canal sealers, a topic of great interest in endodontics. Recently, bioceramic root canal sealers, developed from a modified composition of MTA, have been widely marketed around the world. Against this backdrop, we focused on "Bioactive Glass"<sup>1,2)</sup> which shows similar characteristics to MTA and developed a next-generation bioceramic root canal sealer.

**Characteristics of Bioactive Glass**

Bioactive Glass contains the glass type of Na<sub>2</sub>O-CaO-SiO<sub>2</sub>-P<sub>2</sub>O<sub>5</sub>, as a component of Silica (SiO<sub>2</sub>)

is ≤50mol%. Bioactive Glass has been applied clinically in orthopedic surgery for several decades. When Bioactive Glass is implanted in a defect area close to bone, calcium or a few silicate ions were released from the surface of Bioactive Glass to body fluid, followed by the formation of silica-rich gel on its surface. Silica-rich gel reacts with ions of body fluid, and as a result, hydroxyapatite (HAp)-like crystals form on the surface of Bioactive Glasses (Figure 1). Furthermore, osteoblasts produce new bone in silica-rich gel and Bioactive Glass finally bonds to bone through the interaction between the HAp-like crystals and new bone.





These characteristics show that Bioactive Glass is suitable as a biomaterial for the application to the hard tissues such as dentin or cementum that are similar to bone. It also shows excellent



biocompatibility with soft tissues. Bioactive Glass would contribute to the progress of dental biomaterial used in Endodontics that targets the periapical tissue including hard and soft tissues.

**New Bioceramic root canal sealer**

NISHIKA CANAL SEALER BG was born through the development of biomaterials consisting of Bioactive Glass. "BG" is an abbreviation for "Bioactive Glass". The following are features of NISHIKA CANAL SEALER BG.

**Figure 1****a: Releasing various ions from the surface layer of Bioactive Glass and pH changes****b: HAp formation**

Source: Washio A. et al., Dental Diamond Vol.42-11 No.621 178-183 AUG.2017. (Abridged)

# Features of NISHIKA CANAL SEALER BG

## Physicochemical properties

NISHIKA CANAL SEALER BG is compliant with ISO standard of root canal sealing materials (ISO 6876) (Table 1) and has the property that allows for slight expansion during setting (Figure 2).

After NISHIKA CANAL SEALER BG was immersed in simulated body fluid (SBF), the pH stabilized at approximately 10 and typical spherules of petal-like HAp crystals were observed on the surfaces<sup>3)</sup> (Figure 1). This feature offers excellent sealing ability and biocompatibility.

## Excellent biocompatibility

According to *in vitro* and *in vivo* studies regarding biocompatibility, the excellent

biocompatibility of NISHIKA CANAL SEALER BG was demonstrated. As the important process of wound healing, effects of NISHIKA CANAL SEALER BG on cell migration ability and viability of human periodontal ligament cells (HPDLC) and osteoblast-like cells were examined. As a result, cell migration ability and viability of HPDLC and osteoblast-like cells under NISHIKA CANAL SEALER BG showed no significant difference in comparison with a control<sup>4)</sup>. It is also indicated that HPDLC and osteoblast-like cells migrated and grew in closely contact with the hardened block of NISHIKA CANAL SEALER BG (Figure 3). Furthermore, *in vivo* results demonstrated that NISHIKA CANAL SEALER BG did not inhibit

wound healing of periapical tissue<sup>5)</sup>, showing that NISHIKA CANAL SEALER BG has excellent biocompatibility for periapical tissue.

## Excellent sealing ability

When NISHIKA CANAL SEALER BG was filled in the root canal of human extracted tooth after the irrigation and drying at the same level with clinical procedures, the formation of HAp crystals in the interface between NISHIKA CANAL SEALER BG and root canal wall was observed. Furthermore, the formation of HAp crystals into dentinal tubules and a tag-like structure at the entrance of the tubules were observed. This interaction of NISHIKA CANAL SEALER BG with dentin through the

formation of HAp layer brings its sealing ability<sup>6)</sup> (Figure 4). In a sealing ability test, dye leakage of the root canal filled with gutta-percha point and NISHIKA CANAL SEALER BG by lateral condensation technique was approximately half in comparison with conventional root canal sealers, and the leakage gradually decreased over time<sup>6)</sup>. When a root canal was filled with NISHIKA CANAL SEALER BG by the single-cone technique, the leakage was less than that of the lateral condensation technique. These results showed the excellent sealing ability of NISHIKA CANAL SEALER BG.

**Table 1**  
Physical properties tests (ISO 6876:2012 Requirement)


Flow	28.7mm	Solubility	0.5%
Working time	15 min	Disintegration	None
Setting time	180min	Radiopacity	5mm Al.
Film thickness	27.9μm		

Source: Washio A. et al. The Japanese Journal of Conservative Dentistry, 60(1): 14-21, 2017. (Abridged)



Source: Washio A. et al., Dental Diamond Vol.42-11 No.621 178-183 AUG.2017. (Abridged).


**Figure 3**



HPDL cells and osteoblast-like cells grew around NISHIKA CANALSEALER BG. ( $\times 40$  Magnification)

Source: Washio A. et al., Dental Diamond Vol.42-11 No.621 178-183 AUG.2017. (Abridged).

**Figure 4**



Field emission scanning electron microscope (FE-SEM) photographs of sealer-dentin interfaces of samples (S: Sealer, D: Dentin). Arrows: the formation of tag-like structures in dentinal tubules, arrow heads: hydroxyapatite-like crystals in dentinal tubules.

Source: Washio A. et al., Dental Diamond Vol.42-11 No.621 178-183 AUG.2017. (Abridged).

# Features of NISHIKA CANAL SEALER BG

## Handleability

NISHIKA CANAL SEALER BG is a 2-paste and double syringe model.

By pushing the plunger of a double syringe, it is very easy to dispense 2-paste at a 1:1 ratio. The dispensed paste can be mixed easily and quickly.

Anyone can prepare the same ideal consistency at any time (Figure 5).

## Removability

For the successful retreatment, the removal of root canal sealers from root canal is also one of important properties. *In vitro* study with extracted tooth demonstrated that NISHIKA CANAL SEALER BG was removed completely by the re-preparation and the irrigation, and dentinal tubules of root canal dentin were opened<sup>7)</sup> (Figure


6), showing that NISHIKA CANAL SEALER BG does not inhibit retreatment.

## Precaution

Paste of NISHIKA CANAL SEALER BG tends to get hardened due to heat or moisture absorption. Therefore, NISHIKA CANAL SEALER BG is provided with an aluminum foil resealable bag. It is

recommended that the syringe should be stored in the aluminum foil resealable bag and then stored in a cold place (1-10°C) without freezing.

Regarding mixing paste, stainless-steel spatula may be worn out by its ingredient. Plastic spatula is recommended in order to avoid contaminations.



## References

- 1) Kokubo T. et al.: How useful is SBF in predicting in vivo bone bioactivity? Biomaterial, (27): 2907-2915, 2006.
- 2) Lutz-Christian G. et al.: Bioactive Glass and Glass-Ceramic Scaffolds for Bone Tissue Engineering. Materials, (3): 3867-3910, 2010.
- 3) Washio A. et al.: Physicochemical properties of newly developed bioactive glass cement and its effects on various cells. J Biomed Mater Res B Appl Biomater, 103(2): 373-380, 2015.
- 4) Washio A. et al.: Effects of bioactive glass based sealer on cell migration ability and viability of periodontal ligament cells and osteoblast-like cells. The Japanese Journal of Conservative Dentistry, 60(2), 96-104, 2017.
- 5) Morotomi T. et al.: Effect of newly-developed bioactive glass root canal sealer on periapical tissue of rat's molar. The Japanese Journal of Conservative Dentistry, 60(3): 120 -127, 2017.
- 6) Yoshii S. et al.: Root canal sealing ability of bioactive glass-based sealer and its effects on dentin. The Japanese Journal of Conservative Dentistry, 59(6): 463-471, 2016.
- 7) Washio A. et al.: Evaluation of removability of root canals filled using bioactive glass based sealer. The Japanese Journal of Conservative Dentistry, 60(1): 14-21, 2017.

## Presentation

### Packaging

1 Double syringe [Paste A: 1.5g . Paste B: 1.5g ]

### Composition

Paste A: Fatty acid, Bismuth subcarbonate, Silicon dioxide

Paste B: Magnesium oxide, Purified water,  
Calcium silicate glass (a type of Bioactive Glass), Silicon dioxide, etc.

